


Set up your *Concha Sol* Hearing aids



Download the app

Go to the Apple App Store to download our app “Concha Labs”, using the QR code.



Sign in to activate your hearing aids


Be in a quiet environment, have your iPhone ready, and turn your Bluetooth on!

Hear clearly, live freely


If you would like to view the instructions for use online, please visit conchalabs.com

Concha tip: Follow the app onboarding instructions with care. Please do not put batteries in the hearing aids until instructed.

What's in your box



a.



b.



c.



d.



e.



f.

a. 1 pair of hearing aids with 7mm ear tips in a case

b. Five 6-count packs of disposable zinc-air 312 batteries

c. 1 pair of ear plugs

d. 1 additional pair of 7mm ear tips

e. 2 pairs of 10mm ear tips





f. 1 pack of wax filters

Need help?

 conchalabs.com

- conchalabs.com/conchasol-ifu
- care@conchalabs.com
- +1 (510) 544-6778

Connect with us

-  conchalabs.com
-  [conchalabs](https://twitter.com/conchalabs)
-  [conchalabs](https://www.instagram.com/conchalabs)
-  [Concha Labs](https://www.facebook.com/ConchaLabs)